

NGS analysis for gene regulation and epigenomics

Jacques Serizay, Cyril Matthey-Doret

January 11, 2020

Day 1

Introduction to Hi-C

Regulatory function of chromatin

- Chromatin made up of many DNA-bound proteins
- Dynamic interactions
- Chromatin state regulates gene expression


A picture of the regulatory landscape

- Gene expression tells us what genes are affected during e.g. differentiation, disease...
- Epigenomics can tell us how / why they are activated


Types of chromatin information

We will look at 3 scales of chromatin organization during this workshop:


Adapted from [Illumina](#)



Chromatin conformation capture
(3-C, Hi-C and Capture-C)




Assay for transposase accessible
chromatin (ATAC-Seq)



ChIP-Seq of methylated histones
(Histone methylation)


Capturing chromosome conformation

- 3C methods capture the 3D structure of the genome
- Originally relied on qPCR, Hi-C is the NGS variant




Hi-C experimental process

Interpreting chromosome contact maps



Contact frequency reflects spatial organization

Interpreting chromosome contact maps



- Diagonal gradient due to polymer behavior
- Various patterns correspond to 3D structures

Derivative methods

3C-based methods (digestion + religation):

- MicroC: DNase instead of restriction enzyme
- Promoter capture: Hi-C with capture probes on promoters

Others:

- GAM: Cryoslicing of nuclei followed by sequencing
- SPRITE: Serial dilution and barcoding

From reads to contact maps


Many “end-to-end” pipelines available:

- **nf-core/hic**
- Hi-C pro
- FanC
- hicstuff
- Juicer
- ...

Important to understand individual steps (biases, custom analyses, ...)


Hi-C processing: Overview

General steps common to all Hi-C pipelines




Hi-C processing: Read alignment

- Standard short read mapping
- Separate single-end alignment for forward and reverse
- Longer reads relative to fragment length: More chimeric reads



Hi-C processing: Read alignment


- Solutions generally built in pipelines
 - iterative alignment
 - in-silico fastq digestion
- Mostly important with longer reads



```
truncate read
while not mapped:
    extend read
    mapped = align read
```


Hi-C processing: Merging BAMs

- Merge forward and reverse reads by name
- Filter out ambiguous alignment (low mapQ)
- 1 pair of reads = 1 contact
- The [pairs](#) format is commonly used




Hi-C processing: Fragment attribution

- Assign contacts to discrete segments (bins)
- Usually regular intervals of e.g. 10kb



Hi-C processing: The Matrix

- Sum contact events for each combination of bins
- Most combinations have 0 contacts: “sparse data”
- Storing NxN matrix with mostly 0 is impractical + Use sparse format instead



	0	1	2	3
0				
1		4		
2				
3			2	


0	0	0	0
0	4	0	0
0	0	0	0
0	0	2	0

x	y	v
1	1	4
3	2	2

Hi-C processing: Additional filters

Hi-C specific filters to remove uninformative events

- Self-religating fragments
- Undigested fragments
- Duplicates filter




Hi-C processing: Balancing

Hi-C is susceptible to many biases affecting local coverage such as:

- GC content
- Chromatin accessibility
- Restriction site density

Matrix balancing is a normalization to reduce the impact of those biases.



Hi-C processing: Balancing

Assumption: All bins (regions) have the same contact probability.

- Divide each pixel by row x col iteratively.
- Predefined number of iterations or until convergence.

```
# Simplified ICE procedure
m, n = mat.shape
for _ in range(100):
    for x in range(m):
        row_mean = mat[x, :].mean()
        for y in range(n):
            col_mean = mat[:, y].mean()
            mat[x, y] = mat[x, y] / (row_mean * col_mean)
```

- Generate your own contact map from the reads, normalize and visualize it

Exercises